

Analysis of Dabsyl Amino Acids by High-performance Liquid Chromatography

Introduction

The postcolumn derivatization with HPLC for amino acids analysis has prevailed as a simultaneous separation method of multicomponents by using a selective, high-sensitivity detector and a reversed phase column.

Many results of amino acids analysis using Dabsyl chloride (4-(Dimethylamino)azobenzene-4'-sulfonyl chloride:DABS-Cl) as a reagent for precolumn derivatization have been reported because the amino acids derivatized by Dabsyl chloride are comparatively more stable than by other derivatization agents and for detection, simply a UV/VIS detector can be used with high sensitivity. Here the standard mixture of amino acids was measured by the system using Dabsyl chloride.

Keyword: HPLC, Amino acids, DABS-Cl, DAB Label, 5.0 μm , Dabsylpak II, UV/VIS detector

Experimental

Equipment

Pump: PU-2080

Degasser: DG-2080-53

Gradient unit: LG-2080-02

Autosampler: AS-2055

Column oven: CO-2060

Detector: UV-2070

Conditions

Column: Dabsylpak II (4.6 mmID x 150 mmL, 5 μm)

Pre-column: Dabsylpak II-P (4.6 mmID x 35 mmL, 5 μm)

Eluent A: 20 mM Sodium acetate (pH6.0)

Eluent B: Acetonitrile

Gradient condition: (A/B), 0 min (78/22) \rightarrow 3 min (78/22) \rightarrow

25 min (70/30) \rightarrow 40 min (40/60) \rightarrow

40.1 min (20/80) \rightarrow 45 min (20/80) \rightarrow

45.1 min (78/22) 1 cycle: 55 min

Flow rate: 1.0 mL/min

Column temp.: 45 $^{\circ}\text{C}$

Wavelength: 465 nm

Injection volume: 20 μL

Standard sample: Amino acids mixture 40 pmol each, Type H

Result

Fig. 1 shows the chromatogram of standard mixture of amino acids. 17 kinds of amino acids were separated in 45 minutes.




Fig. 1. Chromatogram of standard mixture of amino acids

1: Aspartic acid, 2: Glutamic acid, 3: Serine, 4: Threonine, 5: Arginine, 6: Glycine, 7: Alanine 8: Proline, 9: Valine, 10: Methionine, 11: Isoleucine, 12: Leucine, 13: Phenylalanine, 14: Cystine 15: Lysine, 16: Histidine, 17: Tyrosine

copyright©JASCO Corporation

Advantages of Dabsylation

- By using the kit (DAB Label), the procedure of intricate preparation of derivatization can be avoided.
- Amino acid composition of 0.5 μg of proteins or peptides can be analyzed with good accuracy and reproducibility.
- 17 kinds of the derivatized amino acids can be separated within 45 minutes.
- The derivatized amino acids are very stable. (for one month at room temperature)
- The derivatized amino acids can be detected by visible light at 465 nm, which allows high-sensitivity analysis without the interference from other components having UV absorption.

DABS-Cl reacts with α-amino group, ε-amino group, phenolic hydroxyl group and imidazole group.

Fig. 2 shows a reaction formula with α-amino group




Fig. 2. Reaction formula with α-amino group