



Recent years have seen explosive growth in the use of organic materials, composites, or organics bound to an inorganic substrate, to solve engineering problems that traditionally were approached only with inorganics. The remarkable speed with which these new technologies are appearing is contrasted by the remarkable scarcity of non-destructive, *in-vitro*, and *in-vivo* techniques for characterizing the biophysical, chemical, and mechanical properties of these often delicate and nano-structured materials. Near-field scanning optical microscopy (NSOM) for quantitative evaluation of surfaces is a practical solution for this type of evaluation.

**Introduction :** Traditional spectroscopy methods utilize 'far-field' light easily propagated through long distances. However, far field light is subject to diffraction effects, whether mirrors or lenses are used. Diffraction effects limit the amount of 'focus' that can be achieved, thus the spatial resolution for traditional microspectroscopy systems is limited to the wavelengths used; i.e. 10 microns for IR, ~1 micron for visible wavelengths.




**Figure 1.** The near field probe.

Near-field scanning optical microscopy (NSOM) is a type of microscopy where a light source with sub-wavelength spatial resolution is used as a scanning probe. The probe is scanned a few nanometers above the sample surface, the piezo-electric stage manipulators maintaining the probe above the surface using a feedback circuit from the probe electronics. A small aperture on the end of a tapered and gold-coated optical fiber is used as the near-field probe (Fig. 1). By illuminating a sample with the "near-field" created by the probe aperture, optical images can be constructed with a spatial resolution well beyond the usual "diffraction limit", typically about 50 nm or greater.



**Figure 2.** Photoluminescence spectra of nanotubes

**Results and Discussion:** Luminescence spectra of nanotubes on a silica surface were collected using an NFS-310 near-field spectrometer. A near-field probe with an aperture of 50 nm was used to illuminate the sample with a 532 nm laser coupled to an optical fiber. An area of 15 x 20 microns was mapped to produce the spectra in Figure 2. Integration of the photoluminescence spectra was used to produce the intensity map in Figure 3. The luminescence intensities associated with discrete nanotubes corresponding to the 'peaks' displayed in the intensity map (Fig. 3).



**Figure 3.** Luminescence intensity map of nanotubes.

**Conclusions:** Near-field scanning optical microscopy offers the ability to produce spectral data with sub-wavelength spatial resolution, beyond the diffraction limits of traditional spectroscopic microscopy methods.

**For More Information Contact:**

Jasco Inc  
1-800-333-5272  
www.jascoinc.com